

# Safety Feature Evaluation Form

## LeEject Syringe

Date: \_\_\_\_\_ Department: \_\_\_\_\_ Occupation: \_\_\_\_\_

Product: \_\_\_\_\_ Number of times used: \_\_\_\_\_

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

### DURING USE:

Agree ..... Disagree

1. The safety feature can be activated using a one-handed technique..... 1 2 3 4 5 N/A
2. The safety feature **does not** obstruct vision of the tip of the sharp..... 1 2 3 4 5 N/A
3. Use of this product requires you to use the safety feature..... 1 2 3 4 5 N/A
4. This product does not require more time to use than a non-safety device..... 1 2 3 4 5 N/A
5. The safety feature works well with a wide variety of hand sizes..... 1 2 3 4 5 N/A
6. The device is easy to handle while wearing gloves..... 1 2 3 4 5 N/A
7. This device **does not** interfere with uses that do not require a needle..... 1 2 3 4 5 N/A
8. This device offers a good view of any aspirated fluid..... 1 2 3 4 5 N/A
9. This device will work with all required syringe and needle sizes..... 1 2 3 4 5 N/A
10. This device provides a better alternative to traditional recapping..... 1 2 3 4 5 N/A

### AFTER USE:

11. There is a clear and unmistakable change (audible or visible) that occurs when the safety feature is activated..... 1 2 3 4 5 N/A
12. The safety feature operates reliably..... 1 2 3 4 5 N/A
13. The exposed sharp is permanently blunted or covered after use and prior to disposal..... 1 2 3 4 5 N/A
14. This device is no more difficult to process after use than non-safety devices..... 1 2 3 4 5 N/A

### TRAINING:

15. The user **does not** need extensive training for correct operation..... 1 2 3 4 5 N/A
16. The design of the device suggests proper use..... 1 2 3 4 5 N/A
17. It is **not** easy to skip a crucial step in proper use of the device..... 1 2 3 4 5 N/A

Of the above questions, which three are the most important to **your** safety when using this product?

Are there other questions which you feel should be asked regarding the safety/utility of this product?



© June 1993, revised August 1998  
Training for Development of Innovative Control Technology Project

# LeEject<sup>®</sup> Syringe & Needle

Dental Safety Device

Distributor Manufacturer:

**Advanced Technology & Capital, Inc.**

5 Sylvan Avenue, Englewood Cliffs, NJ 07632

201-944-8400

info@LeEject.com